

130mg Clean Screen Xcel I Column
 Part #: CSXCE106 6 mL - 130 mg Cartridge
 Part #: CSXCE103 3 mL - 130 mg Cartridge



1. SAMPLE PREPARATION

(Hydrolysis Step)

To 1-5 mL urine sample add 1-2 mL of 0.1M acetate buffer (pH= 5.0) containing 5,000 units/mL β -glucuronidase.

Optionally, add 500 μ L of acetate buffer and 25 μ L of concentrated β -glucuronidase.

Add appropriate volume and concentration internal standards.

Vortex and heat for 1-2 hours at 65 °C.

Allow sample to cool.

Do not adjust pH- sample is ready to be added to extraction column.

2. APPLYING SAMPLE TO COLUMN

Load sample directly to column without any preconditioning.

Pull sample through at a rate of 1-2 mL/ minute.

Dry column thoroughly under vacuum (10 mm Hg) or positive pressure (~ 80-100 psi) for 1 minute.

3. WASH

Wash sample with 1 mL of methylene chloride.

Dry column thoroughly under vacuum (10 mm Hg) or positive pressure (~ 80-100 psi) for a minimum of 5-10 minutes.



NOTE 2: (It is important to dry the column thoroughly to achieve the highest recovery of all compounds. Any residual moisture will slow down the drying of the elution solvents prior to derivatization for GC/MS analysis. Also, any residual moisture could reduce the reactivity of the derivatization agent resulting in low GC/MS sensitivity.)

4. ELUTION

Elute samples with 1 mL ethyl acetate/ ammonium hydroxide (98/2) Evaporate fraction to complete dryness under stream of dry air or nitrogen at ~ 35 °C.

GC/MS Analysis

It is recommended to add 50 µL of Ethyl Acetate to 50 µL of derivatization agent and react at ~70 °C for 15 minutes. Inject 1-2 µL of cooled (50:50) solution in the GC/MS system for analysis.

LC/MS Analysis

Reconstitute in methanol or appropriate mobile phase.

Benzodiazepine Analytes Extracted

Diazepam	Clonazepam	Alprazolam	Midazolam
Nordiazepam	7-aminoclonazepam	α-OH Alprazolam	α-OH Midazolam
Temazepam	Oxazepam	Lorazepam	Flurazepam