



**COCAINE, BENZOYLECGONINE, AND COCAETHYLENE IN BLOOD,
PLASMA/SERUM, URINE AND TISSUE USING:
200 mg CLEAN SCREEN® EXTRACTION COLUMN**
Part #: ZSDAU020
LC-MSMS

1. PREPARE SAMPLE:

To 1 mL of 100 mM phosphate buffer (pH 6) add internal standard.*

Add 1 mL of whole blood, serum/ plasma, urine, or 1g tissue homogenate (1:4). Add 2 mL of 100 mM phosphate buffer (pH 6).

Vortex and centrifuge as appropriate.

2. CONDITION COLUMN:

1 x 3 mL CH₃OH

1 x 3 mL D.I. H₂O

1 x 1 mL 100 mM phosphate buffer (pH 6).

Note: aspirate at < 3 inches Hg to prevent sorbent drying out

3. APPLY SAMPLE:

Load sample at 1-2 mL / minute.

4. WASH COLUMN:

1 x 3 mL DI H₂O

1 x 3 mL 100 mM hydrochloric acid.

1 x 3 mL CH₃OH

Dry column (5 minutes at > 10 inches Hg).

5. ELUTE COCAINE, BENZOYLECGONINE, COCAETHYLENE:

1 x 3 mL ethyl acetate; acetonitrile: ammonia (78: 20: 2 v/v)

Or

1 x 3 mL CH₂Cl₂/ IPA/ ammonia (78:20:2 v/v)

Collect eluate at 1-2 mL /minute.

6. EVAPORATION:

Evaporate eluates to dryness under a gentle stream of nitrogen.

7. RECONSTITUTE sample in 50 µL of CH₃OH.

Inject 5 µL.

7. INSTRUMENT CONDITIONS:

Column: 50 x 2.1 mm (3 µm) Selectra[®] Phenyl (UCT, Inc.,)

MOBILE PHASE:

<u>Time</u>	<u>% Acetonitrile</u>	<u>%0.1% Formic acid</u>
0	33	67
5	33	67

Flowrate: 0.35 mL/ minute

Column Temperature: ambient

Detector: API 2000 MS/MS


<u>Compound</u>	<u>MRM Transition</u>	<u>Cerilliant #</u>
Cocaine	304.2/ 182.3	C-008
* Cocaine-D3	307.2/ 185.2	C-009
Benzoylecgonine	290.1/168.0	B-004
*Benzoylecgonine-D8	298.2/171.3	B-013
Cocaethylene	318.2/ 196.2	C-010
*Cocaethylene-D8	326.2/204.2	C-024

CHROMATOGRAM SHOWING:

Cocaine/ Cocaine-D3

Benzoyllecginine/ Benzoyllecgonine-D8

Cocaethylene/ Cocaethylene-D8



Recovery: > 90% (n=10)

LOD: 10 ng/ mL