



## Analysis of 1,4-Dioxane in 500 mL of Drinking Water by Gas Chromatography-Mass Spectrometry SIM Detection

Part Number: EU52112M6, 2 grams activated carbon in a 6 mL cartridge

1,4-dioxane (CAS 123-91-1), often referred to as dioxane, is a highly soluble non-biodegradable ether. This compound is used as solvent stabilizer to prevent the breakdown of chlorinated solvents.

### 1) Sample Preparation:

- a) To dechlorinate the water sample add 25 mg of sodium sulfite, **then** acidify with 0.5 g of sodium bisulfate and swirl until dissolved.

### 2) Condition Cartridge:

- a) Place a EU52112M6 cartridge on a single or multi-station vacuum manifold or automated extraction system
- b) Add 3 mL of dichloromethane (methylene chloride) to the cartridge and draw through
- c) Add 3 mL of methanol and draw through
- d) Repeat step c)
- e) Add 3 mL of DI water and draw through
- f) Repeat step e) 4 times

### 3) Load Sample:

- a) Add Internal standard \* to the sample
- b) Add sample to the cartridge. Adjust for a flow rate of 5 - 10 mL/minute
- c) Dry cartridge by drawing air at full vacuum for 10 minutes

### 4) Elution:

- a) Place a clean collection vial in the vacuum manifold
- b) Add 3 mL of dichloromethane to the cartridge. Slowly draw through
- c) Add a second 6 mL aliquot of dichloromethane and draw through
- d) Add surrogate\*\* (0.5 ppm) to the extract and bring to 10 mL final volume with dichloromethane
- e) Dry extract by adding 2 g anhydrous sodium sulfate ( $\text{Na}_2\text{SO}_4$ ) and vortex to mix the extract

**5) Quantitate:**

a) GC-MS-SIM (electron ionization)

**6) Instrument & Conditions:**

**Column:** Varian CP-Select 624 CB (6% cyanopropyl phenyl, 94% PDMS, 30 m 0.25mm x 1.4 um column) or equivalent

**Injector:** 200°C (Splitless mode)

**Injector Volume:** 1 µL

**Flow:** 1 mL minute

**Oven:** 30°C for 1 minute, 90°C at 7°C/minute, 200°C at 20°C/minute for 3 minutes

**MS:** quadrupole MS

**SIM MODE**

**Dwell time 100 µs**

**Emission current 100 µA**

**\*\*\* Quantitation ions**

<b>Segment 1 (0-8 minutes)</b>	<b>**THF-d<sub>8</sub></b>	<b>m/z 40***, 78, 80</b>
<b>Segment 2 (8-16.5 minutes)</b>	<b>1, 4 dioxane</b>	<b>m/z 58, 88***</b>
	<b>*1, 4 dioxane-d<sub>8</sub></b>	<b>m/z 62, 64, 96***</b>

R 1.0 20090114

\*Grinnett, Paul E., Munch, Jean W., Method Development for the Analysis of 1,4-Dioxane in Drinking Water Using Solid-Phase Extraction and Gas Chromatography-Mass Spectrometry, Journal of Chromatographic Science, Vol. 47, January 2009. Office of Research and Development, National Exposure Laboratory, Cincinnati, OH.

DCN-904220-135