



**EXTRACTION OF ANTICOAGULANTS FROM URINE  
USING CEC30 SORBENT IN 96 WELL PLATE FORMAT  
LC-MS/MS**

April 13, 2011

**1. PREPARE SAMPLE:**

To 200  $\mu$ L of 100 mM phosphate buffer (pH 6) add internal standard.\*

Vortex mix

Add 100-500  $\mu$ L of urine and vortex mix

Add 200  $\mu$ L of 100 mM phosphate buffer and vortex mix.

Centrifuge as appropriate

**2. CONDITION 96 WELL PLATES:**

1 x 1 mL  $\text{CH}_3\text{OH}$

1 x 1 mL D.I.  $\text{H}_2\text{O}$

1 x 0.5 mL 100 mM phosphate buffer (pH 6)

**Note:** aspirate at < 3 inches Hg to prevent sorbent drying out

**3. APPLY SAMPLE:**

Load sample at 1-2 mL / minute.

**4. WASH 96 WELL PLATES:**

1 x 1 mL 100 mM phosphate buffer (pH 6).

Dry column (5 minutes at > 10 inches Hg).

1 x 1 mL Hexanes

Dry column (5 minutes at > 10 inches Hg).

**5. ELUTE ANTICOAGULANTS:**

2 x 2 mL of ethyl acetate: methanol (75:25)

Collect eluate at 1-2 mL /minute.

## 6. EVAPORATION:

Evaporate eluates to dryness under a gentle stream of nitrogen (<40°C)

### **Alternative:**

Transfer eluate to a clean glass tube with the aid of a disposable glass pipette.

Wash 96 well collection rack with 1 x 1 mL of elution solvent.

Transfer washing to glass tube containing eluate

Evaporate eluates to dryness under a gentle stream of nitrogen (<40°C)

## 7. LC-MS/MS

Reconstitute sample in 100 µL CH<sub>3</sub>OH Inject 5 µL.

**Column:** 50 x 2.1 (3 µm) SELECTRA Phenyl (UCT, Inc.)

### MOBILE PHASE

<u>Time / min</u>	<u>% Acetonitrile</u>	<u>%0.1% Formic Acid</u>
0	60	40
10	60	40


**Flowrate:** 0.3 mL/minute.

**Column Temperature:** 40°C.

**Detector:** API 4000 QTRAP MS/MS

<u>Compound</u>	<u>(+) MRM Transition</u>	<u>(-) MRM Transition</u>
Warfarin	309.1/163	306.9/161
*p-Chlorowarfarin	343.1/163	340.9/160.0
Broadificoum	524.3/256	520.8/78
Bromodione	527.8/250.5	524.7/249
Coumatetrayl	293.8/175.0	
Flocoumafen	542.4/159.1	

# CHROMATOGRAM OF ANTICOAGULANTS (500 ng) extracted from urine (500 $\mu$ L)



Recovery > 90%

Matrix effects < 10%

DCN-113140-215